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03.07.2017

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Report Form Manufacturer's Trend Report

Medical Devices Vigilance System (MEDDEV 2.12/1 rev 7)

v.12/11

1. Administration Information	
Recipient (Name of National Competent Authority NCA) Danish Medicines Agency	
Address of National Competent Authority Axel Heides Gade 1, DK-2300 Copenhagen	
Date of this report 2017-07-03	
Reference number assigned by the manufacturer 2016050310	
Reference number assigned by NCA 2016050310	
Type of report <input type="checkbox"/> Trend Initial <input type="checkbox"/> Trend Follow up <input checked="" type="checkbox"/> Trend Final	
Do these incidents / trend represent a serious public health threat? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	
Identify to what other NCAs this report was also sent N/A	
2. Information on submitter of the report	
Status of submitter <input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Authorised Representative within EEA, Switzerland and Turkey <input type="checkbox"/> Others: (identify the role) :	
3. Manufacturer information	
Name Novo Nordisk A/S	
Contact name Erik Steffensen	
Address Krogshøjvej 55	
Postcode DK-2880	City Bagsværd
Phone +45 3079 3019	Fax N/A
E-mail esff@novonordisk.com	Country Denmark
4. Authorised Representative information	

Name N/A	
Contact name	
Address	
Postcode	City
Phone	Fax
E-mail	Country
5. Submitter's information (if different from section 3 or 4)	
Submitter's name N/A	
Contact name	
Address	
Postcode	City
Phone	Fax
E-mail	Country
6. Medical Device Information	
Class <input type="checkbox"/> AIMD Active Implants <input type="checkbox"/> MDD Class III <input checked="" type="checkbox"/> MDD Class IIb <input type="checkbox"/> MDD Class IIa <input type="checkbox"/> MDD Class I <input type="checkbox"/> IVD Annex II List A <input type="checkbox"/> IVD Annex II List B <input type="checkbox"/> IVD Devices for self-testing <input type="checkbox"/> IVD General	
Nomenclature system (preferable GMDN) GMDN	Nomenclature code 45771
Nomenclature text Basic insulin pen, replaceable cartridge	
Commercial name/ brand name / make NovoPen [®] Echo [®] and NovoPen [®] 5	
Model number(s) or Family name N/A	Catalogue number(s) N/A
Serial number range (if applicable) N/A	Lot/batch number range(if applicable) N/A
Software version number (if applicable) N/A	
Accessories / associated devices (if applicable) N/A	

Notified Body (NB) ID – Number 0088
7. Information on Trend Report
Date the trend was identified 2017-06-09
Description narrative for identified trend Increased complaint trend related to damaged cartridge holders on NovoPen® Echo® and NovoPen® 5
Time period of trend analysis 2016-2017
Established trigger level Statistically significant increase in complaints related to damaged cartridge holder
Have any of the trended events been submitted individually as reportable events under vigilance? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
If yes, please list how many and to which Competent Authority From data extracted on the 2 nd June 2017, three (3) cases were identified as reported to competent authorities. Two (2) of the cases have been reported to the United Kingdom (Medicines & Healthcare products Regulatory Agency; MHRA) and one (1) has been reported to France (L'Agence nationale de sécurité du médicament et des produits de santé; ANSM).
8. Manufacturer's preliminary comments
Manufacturer's preliminary analysis into causes of trend See investigation under section 9.
Initial corrective actions / preventive actions implemented by the manufacturer Production was changed back to the original cartridge holder material, implemented 2016-09-01.
Expected date of next report N/A
9. Results of manufacturer's final investigation into trend
The manufacturer's trend analysis results Novo Nordisk has evaluated that a NovoPen® Echo® or a NovoPen® 5 device having a damaged cartridge holder could in some cases result in the patient experiencing high blood glucose. Back in July 2014, a re-designed cartridge holder for NovoPen® Echo® and NovoPen® 5 was implemented to improve the robustness towards the increased drop-test requirements in the updated 2012 version of the ISO11608-1 standard for pen-injectors. Unfortunately, the re-designed cartridge holder can become damaged during use under special circumstances, during which the snaps keeping the cartridge holder in place will crack and/or brake off (see figure 1, next page). Customer complaints were received in the beginning of 2016 and the increased complaint trend was reported to DKMA April 2016 (Ref. 2016050310). As the re-designed cartridge holder in combination with inappropriate handling could lead to a defective product, Novo Nordisk in April 2016 decided to change production back to the original cartridge holder so that going forward all products would be produced with the original cartridge holder. This was implemented 1 st of September 2016. The <i>cause of failure</i> was identified to be exposure of the device to certain chemicals and oils during use in combination with the changed cartridge holder design (new material). In the instruction for use the following is described in relation to cleaning: <ul style="list-style-type: none"> • Only clean your pen with mild detergent on a moistened cloth. • Do not wash, soak, or lubricate your pen, and do not use products containing bleaching agents, such as chlorine, iodine, or alcohol to clean your pen. This may damage it. The <i>direct cause</i> was identified to be a result of the new selected material in the re-designed cartridge holder (Polycarbonate (PC)/acrylonitrile butadiene styrene (ABS)) and a thinner geometry of the new design, lending itself for induction of Environmental Stress Cracking in combination with unintended handling by the patients. The original material PC/Polybutylene terephthalate (PBT) and a thicker geometry under the snap detail have not caused a similar issue for devices on the market. The <i>root cause</i> was identified to be lack of adequate Design Verification Test.

Evaluation of criticality

Three types of faults can occur: A) 'Snap cracked', B) 'One snap broken off' and C) 'Both snaps broken off' (figure 1).

Figure 1. Three types of faults can occur: 1) 'Snap cracked', 2) 'One snap broken off and 3) 'Both snaps broken off'.

The three types of faults have been evaluated for criticality according to severity classes used by Novo Nordisk (see Table 1).

Non-safety No medical consequence		Safety With medical consequences		
S1	S2	S3	S4	S5
Dissatisfaction of quality expectations	Discomfort or loss of functionality or quality	Non-serious adverse event	Serious adverse event without fatal outcome	Serious adverse event with fatal outcome

Table 1. Severity Classification – definition.

All three different fault types can emerge during handling of the cartridge holder at the snaps due to the chemical resistance of the material when combined with unintended use of cleaning agents and/or household substances (sunscreen, food grease, etc.).

Descriptions of the snap feature cracking are documented in the Hazard and Operability (HAZOP) Analysis of NovoPen® Echo®. The description covers the error caused by the combination of force and chemicals. It should be noted that the faults are likely progressive and develop from a cracked snap to one snap broken off and finally to both snaps broken off.

In table 2 below, the clinical consequence and the estimated severity for the three fault types shown in figure 1 are summarised.

Fault type	Clinical consequence and estimated severity
Snaps on Cartridge Holder are cracked	Hyperglycaemia S3 / Hypoglycaemia S3 or S4
One snap is broken off	Hyperglycaemia S3 to S5 / Hypoglycaemia S3 or S4
Both snaps are broken off	Hyperglycaemia S3

Table 2. Clinical consequences and estimated severity for the three fault types shown in Figure 1.

The estimated probability and severity for the different fault scenarios are shown in table 3 below.

Hazardous situation	Probability and severity
Over-dose from alignment	10^{-4} , S3 / 10^{-5} , S4
Under-dose following over-dose	10^{-5} , S3
Multiple no-doses	10^{-9} , S5
No treatment alternative	10^{-4} , S3

Table 3. The estimated risk that a patient will experience an event over the lifetime of the pen.

Consequently, based on the risk analysis the identified faults vary in severity from S3 to S5 scenarios and can cause hyper- or hypoglycaemia.

NovoPen® Echo® and NovoPen® 5 on the market

In total around 300,000 affected NovoPen® Echo® and 200,000 NovoPen® 5 devices are on the market, see enclosures for a list of affected batches.

Affected countries with NovoPen® Echo®

Australia, Austria, Belgium, Bulgaria, Canada, Croatia, Czech Republic, Denmark, Finland, France, Georgia, Germany, Greece, Hong Kong, Hungary, Iceland, Ireland, Italy, Kyrgyzstan, Netherlands, New Zealand, Norway, Poland, Portugal, Romania, Russia, Serbia, Singapore, Slovakia, Slovenia, South Africa, Spain, Sweden, Switzerland, Taiwan, Ukraine, United Kingdom and United States.

Affected countries with NovoPen® 5

Austria, Canada, Denmark, Finland, Germany, Hong Kong, Ireland, Netherlands, Sweden, Switzerland and United Kingdom.

In April 2016 when the increased complaint trend was identified only a few adverse events were reported and there was no clear correlation between a cracked cartridge holder and adverse events. However, in April-May 2017 a large increase in the number of complaints was seen and it was decided to perform a new analysis of the adverse events and update the Health Hazard Evaluation.

Based on the actual data from the market we know that S3 adverse events occur at a rate below 10^{-3} , i.e. this correlates to a maximum of 500 S3 adverse events if no action is taken (assuming 500,000 pens on the market). We can expect in total less than 0.5 S4 events and 0.005 S5 events (meaning highly unlikely) if the affected devices stay on the market until end of expiry.

Customer complaints

Since the introduction of the re-designed cartridge holder, 972 verified complaints have been received (as of 31st of May 2017, see figure 2), comprising both NovoPen® 5 and NovoPen® Echo® (391 and 581 complaints, respectively). For NovoPen® Echo and NovoPen® 5 the return rate of complaint samples is above 90 %. Thus, the verified faults give a good representation of the total complaint picture.


Figure 2. Development in global number of customer complaints on NovoPen® 5 and NovoPen® Echo® where the cartridge holder on the returned sample is damaged.

Causality assessment of AE cases

End of May 2017, 67 customer complaints - where the returned pen had a re-designed cartridge holder that was damaged - had an associated reported adverse event. The clinical consequence reported in the majority of cases was non-serious hyperglycaemia (62 cases). In the remaining cases serious hyperglycaemia (2 cases), serious hypoglycaemia (1 case), non-serious malaise (1 case) and non-serious cataract (1 case) were reported. The result of the causality assessment is shown in table 4 below.

Causality	Non-serious cases	Serious cases
Possible ¹	18	0
Unknown ²	34	1
Unlikely ³	12	2

Table 4. Causality assessment of AE cases.

1. Possible: The number and/or specificity of the confounding factors suggest a possible relationship.
2. Unknown: The confounding factors are pointing in different directions or too little information was available.
3. Unlikely: The number and/or specificity of the confounding factors suggest that a causal relationship is unlikely.

A total of three serious cases were received. In two of the serious cases it was evaluated that the experienced adverse events (AEs) were not related to the observed broken/cracked snaps on the cartridge holder while in the third limited and conflicting information was given. In the first case the reported AE did not match the medical consequence listed in the result code for the relevant fault and in the second case the patient consistently re-used needles and left them attached to the pen in between injections. Both needle re-use and leaving the needle attached to the pen are confounding factors which is more likely to cause the reported hyperglycaemic event. In the last case conflicting information related to the reported event and the patients' blood glucose measurements made the outcome of the assessment unknown.

In 18 non-serious cases the data available suggested that there was a causal relationship between the reported AEs and the observed fault(s) on the returned pen. In all of the 18 cases the patient experienced a hyperglycaemic event and in 5 of the 18 cases, the case was medical confirmed by a health care professional. Several factors which suggested a causal relationship or ruled out other potential explanations were present in the cases, please note that more than one factor can be present in one case.

NovoPen® Echo® used in clinical trial NN1218-4101

Clinical trial NN1218-4101 is on-going and has used pens from one of the affected batches. NN1218-4101 is a phase 3b clinical trial on efficacy and safety of faster-acting insulin Aspart compared to NovoRapid®, both in combination with insulin degludec in children and adolescents with type 1 diabetes. NovoPen Echo® is used to administer faster-acting insulin Aspart and insulin Aspart.

A total of 8 complaints are concluded as related to one or both snaps broken off on the cartridge holder. None of these are reported as related to adverse events.

The following countries are affected:

Bulgaria, Czech Republic, Estonia, Finland, Germany, India, Israel, Italy, Japan, Latvia, Lithuania, Poland, Russia, Serbia, Turkey, Ukraine and United States.

There are 433 pens in stock at the central depot and a total of 178 pens at local depots and 1399 pens out Novo Nordisk control. All pens come from a single batch.

Health Hazard Evaluation

The HHE is a class II: A situation in which the use of or exposure to a violative product may cause temporary or medically reversible adverse health consequences, or where the probability of serious adverse health consequences is remote.

Remedial action / corrective action / preventive action / Field Safety Corrective Action

Corrective actions have been implemented in production ensuring that all new NovoPen® Echo® and NovoPen® 5 devices have a cartridge holder without this potential issue.

The consequence of a patient level recall will be a stock out situation for at least 4-6 weeks in certain markets which will put patients in a poor situation as it can have health consequences to switch to another injection device or another type of medication. To ensure patient safety it is decided to:

- Provide affected patients with a cartridge holder of the original and more robust material, which they can easily replace themselves. Novo Nordisk is establishing the logistics, enabling courier service distribution of exchange cartridge holders to all affected patients globally. The cartridge holder is a part of the device and the patient is used to detach the cartridge holder, load a new insulin cartridge and attach the cartridge holder again.
- Recall affected devices from wholesalers, pharmacies, hospitals and doctors.

For clinical trial NN1218-4101, investigators and patients in the possession of affected NovoPen® Echo® devices will be contacted to have their devices replaced.

Time scheduled for the implementation of the identified actions

2017-07-05: Initiate Field Safety Corrective Actions outlined above in agreement with local health authorities

Final comments from the manufacturer

N/A

Further investigation

N/A

10. The medical device has been distributed to the following Countries

Within EEA, Switzerland and Turkey:

<input checked="" type="checkbox"/> AT	<input checked="" type="checkbox"/> BE	<input checked="" type="checkbox"/> BG	<input checked="" type="checkbox"/> CH	<input type="checkbox"/> CY	<input checked="" type="checkbox"/> CZ	<input checked="" type="checkbox"/> DE	<input checked="" type="checkbox"/> DK	<input type="checkbox"/> EE	<input checked="" type="checkbox"/> ES
<input checked="" type="checkbox"/> FI	<input checked="" type="checkbox"/> FR	<input checked="" type="checkbox"/> GB	<input checked="" type="checkbox"/> GR	<input checked="" type="checkbox"/> HU	<input checked="" type="checkbox"/> IE	<input checked="" type="checkbox"/> IS	<input checked="" type="checkbox"/> IT	<input type="checkbox"/> LI	<input type="checkbox"/> LT
<input type="checkbox"/> LU	<input type="checkbox"/> LV	<input type="checkbox"/> MT	<input checked="" type="checkbox"/> NL	<input checked="" type="checkbox"/> NO	<input checked="" type="checkbox"/> PL	<input checked="" type="checkbox"/> PT	<input checked="" type="checkbox"/> RO	<input checked="" type="checkbox"/> SE	<input checked="" type="checkbox"/> SI
<input checked="" type="checkbox"/> SK	<input type="checkbox"/> TR								

Candidate Countries:

HR

All EEA, Candidate Countries, Switzerland and Turkey

All countries:

Marketed products:

Australia, Austria, Belgium, Bulgaria, Canada, Croatia, Czech Republic, Denmark, Finland, France, Georgia, Germany, Greece, Hong Kong, Hungary, Iceland, Ireland, Italy, Kyrgyzstan, Netherlands, New Zealand, Norway, Poland, Portugal, Romania, Russia, Serbia, Singapore, Slovakia, Slovenia, South Africa, Spain, Sweden, Switzerland, Taiwan, Ukraine, United Kingdom and United States.

Clinical trial:


Bulgaria, Czech Republic, Estonia, Finland, Germany, India, Israel, Italy, Japan, Latvia, Lithuania, Poland, Russia, Serbia, Turkey, Ukraine and United States.

11. Comments

Submission of this report does not, in itself, represent a conclusion by the manufacturer and / or authorized representative or the National Competent Authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical

device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.

I affirm that the information given above is correct to the best of my knowledge.

.....
Erik Steffensen Bøgsværd 2017-07-03

Report Form

Field Safety Corrective Action

Medical Devices Vigilance System

(MEDDEV 2.12/1 rev 7)

new case, keep base data

Version 2.7en
2012-12-03

1 Administrative information	
To which NCA(s) is this report being sent?	
To the NCAs in the countries where the affected batches are shipped.	
Type of report	
<input checked="" type="radio"/> Initial report <input type="radio"/> Follow-up report <input type="radio"/> Final report	
Date of this report	
2017-07-03	
Reference number assigned by the manufacturer	
2016050310	
FSCA reference number assigned by NCA	
2016050310	
Incidence reference number assigned by NCA	
Name of the co-ordinating NCA Competent Authority (if applicable)	

2 Information on submitter of the report

Status of submitter

- Manufacturer
 Authorised Representative within EEA and Switzerland
 Others: (identify the role)

3 Manufacturer information

new

Name	
Novo Nordisk A/S	
Contact Name	
Erik Steffensen	
Address	
Krogshøjvej 55	
Postcode	City
DK-2880	Bagsværd
Phone	Fax
+45 3079 3019	
E-mail	Country
esff@novonordisk.com	DK - Denmark

4 Authorised Representative Information

new

Name	
Contact Name	
Address	
Postcode	City
Phone	Fax
E-mail	Country DE - Germany

5 National contact point information

new

National contact point name Novo Nordisk A/S	
Name of the contact person Erik Steffensen	
Address Krogshøjvej 55	
Postcode DK-2880	City Bagsværd
Phone +45 3079 3019	Fax N/A
E-mail esff@novonordisk.com	Country DK - Denmark

6 Medical device information

new

Class

- AIMD Active implants
 MDD Class III
 MDD Class IIb
 MDD Class IIa
 MDD Class I
- IVD Annex II List A
 IVD Annex II List B
 IVD Devices for self-testing
 IVD General

Nomenclature system (preferable GMDN)

GMDN

Nomenclature code

45771

Nomenclature text

Basic insulin pen, replaceable cartridge

Commercial name/ brand name / make

NovoPen® Echo® and NovoPen® 5

Model number**Catalogue number****Serial number(s)****Lot/batch number(s)**

See enclosures 1-2.

Device Mfr Date**Expiry date****Notified Body (NB) ID-number**

0088

Accessories / associated devices (if applicable)**Software version number (if applicable)**

7 Description of the FSCA

Background information and reason for the FSCA

Novo Nordisk has evaluated that a NovoPen® Echo® or a NovoPen® 5 device having a damaged cartridge holder could in some cases result in the patient experiencing high blood glucose.

Back in July 2014, a re-designed cartridge holder for NovoPen® Echo® and NovoPen® 5 was implemented to improve the robustness towards the increased drop-test requirements in the updated 2012 version of the ISO11608-1 standard for pen-injectors. Unfortunately, the re-designed cartridge holder can become damaged during use under special circumstances, during which the snaps keeping the cartridge holder in place will crack and/or break off. Customer complaints were received in the beginning of 2016 and the increased complaint trend was reported to DKMA April 2016 (Ref. 2016050310). As the re-designed cartridge holder in combination with inappropriate handling could lead to a defective product, Novo Nordisk in April 2016 decided to change production back to the original cartridge holder so that going forward all products would be produced with the original cartridge holder. This was implemented 1st of September 2016.

The cause of failure was identified to be exposure of the device to certain chemicals and oils during use in combination with the changed cartridge holder design (new material). In the instruction for use the following is described in relation to cleaning:

- Only clean your pen with mild detergent on a moistened cloth.
- Do not wash, soak, or lubricate your pen, and do not use products containing bleaching agents, such as chlorine, iodine, or alcohol to clean your pen. This may damage it.

The direct cause was identified to be a result of the new selected material in the re-designed cartridge holder (Polycarbonate (PC)/acrylonitrile butadiene styrene (ABS)) and a thinner geometry of the new design, lending itself for induction of Environmental Stress Cracking in combination with unintended handling by the patients. The original material PC/Polybutylene terephthalate (PBT) and a thicker geometry under the snap detail have not caused a similar issue for devices on the market.

The root cause was identified to be lack of adequate Design Verification Test.

EVALUATION OF CRITICALITY

Three types of faults can occur with the cartridge holder:

- A) 'Snap cracked',
- B) 'One snap broken off' and
- C) 'Both snaps broken off' (see figure below).

The three types of faults have been evaluated for criticality according to severity classes used by Novo Nordisk.

Non-safety (No medical consequence):

- Dissatisfaction of quality expectations = S1
- Discomfort or loss of functionality or quality = S2

Safety With medical consequences:

- Non-serious adverse event = S3
- Serious adverse event without fatal outcome = S4
- Serious adverse event with fatal outcome = S5

All three different fault types can emerge during handling of the cartridge holder at the snaps due to the chemical resistance of the material when combined with unintended use of cleaning agents and/or household substances (sunscreen, food grease, etc.).

Descriptions of the snap feature cracking are documented in the Hazard and Operability (HAZOP) Analysis of NovoPen® Echo®. The description covers the error caused by the combination of force and chemicals. It should be noted that the faults are likely progressive and develop from a cracked snap to one snap broken off and finally to both snaps broken off.

Below, the clinical consequence and the estimated severity for the three fault types are summarized.

Fault type, Clinical consequence and estimated severity:

- Snaps on Cartridge Holder are cracked: Hyperglycemia S3 / Hypoglycemia S3 or S4
- One snap is broken off: Hyperglycaemia S3 to S5 / Hypoglycaemia S3 or S4
- Both snaps are broken off: Hyperglycaemia S3

The probability and severity for the different fault scenarios are listed below.

Hazardous situation, Probability and severity:

- Over-dose from alignment: 10E-4, S3 / 10E-5, S4
- Under-dose following over-dose: 10E-5, S3
- Multiple no-doses: 10E-9, S5
- No treatment alternative: 10E-4, S3

Consequently, based on the risk analysis the identified faults vary in severity from S3 to S5 scenarios and can cause hyper- or hypoglycaemia.

NovoPen® Echo® and NovoPen® 5 on the market:

In total around 300000 affected NovoPen® Echo® and 200000 NovoPen® 5 devices are on the market, see enclosures for a list of affected batches.

Affected countries with NovoPen® Echo®

Australia, Austria, Belgium, Bulgaria, Canada, Croatia, Czech Republic, Denmark, Finland, France, Georgia, Germany, Greece, Hong Kong, Hungary, Iceland, Ireland, Italy, Kyrgyzstan, Netherlands, New Zealand, Norway, Poland, Portugal, Romania, Russia, Serbia, Singapore, Slovakia, Slovenia, South Africa, Spain, Sweden, Switzerland, Taiwan, Ukraine, United Kingdom and United States.

Affected countries with NovoPen® 5

Austria, Canada, Denmark, Finland, Germany, Hong Kong, Ireland, Netherlands, Sweden, Switzerland and United Kingdom.

In April 2016 when the increased complaint trend was identified only a few adverse events were reported and there was no clear correlation between a cracked cartridge holder and adverse events. However, in April-May 2017 a large increase in the number of complaints was seen and it was decided to perform a new analysis of the adverse events and update the Health Hazard Evaluation.

Based on the actual data from the market we know that S3 adverse events occur at a rate below 10-3, i.e. this correlates to a maximum of 500 S3 adverse events if no action is taken (assuming 500,000 pens on the market). We can expect in total less than 0.5 S4 events and 0.005 S5 events (meaning highly unlikely) if the affected devices stay on the market until end of expiry.

NovoPen® Echo® used in clinical trial NN1218-4101:

Clinical trial NN1218-4101 is on-going and has used pens from the affected batches. NN1218-4101 is a phase 3b clinical trial on efficacy and safety of faster-acting insulin Aspart compared to NovoRapid®, both in combination with insulin degludec in children and adolescents with type 1 diabetes. NovoPen Echo® is used to administer faster-acting insulin Aspart and insulin Aspart.

8 complaints are concluded as related to one or both snaps broken off on the cartridge holder. None of these are reported as related to adverse events.

The following countries are affected:

Bulgaria, Czech Republic, Estonia, Finland, Germany, India, Israel, Italy, Japan, Latvia, Lithuania, Poland, Russia, Serbia, Turkey, Ukraine and United States.

There are 433 pens in stock at the central depot and a total of 178 pens at local depots and 1399 pens out Novo Nordisk control. All pens come from a single batch.

CUSTOMER COMPLAINTS

Since the introduction of the re-designed cartridge holder, 972 verified complaints have been received (as of 31 May 2017), comprising both NovoPen® 5 and NovoPen® Echo® (391 and 581 complaints, respectively). For NovoPen® Echo and NovoPen® 5 the return rate of complaint samples is above 90 %. Thus, the verified faults give a good representation of the total complaint picture.

Causality assessment of AE cases:

End of May 2017, 67 customer complaints - where the returned pen had a re-designed cartridge holder that was damaged - had an associated reported adverse event. The clinical consequence reported in the majority of cases was non-serious hyperglycaemia (62 cases), in the remaining cases serious hyperglycaemia (2 cases), serious hypoglycaemia (1 case), non-serious malaise (1 case) and non-serious cataract (1 case) were reported. The result of the causality assessment is shown below.

Causality: Non-serious cases // Serious cases

- Possible (1): 18 // 0
- Unknown (2): 34 // 1
- Unlikely (3): 12 // 2

1. Possible: The number and/or specificity of the confounding factors suggest a possible relationship.
2. Unknown: The confounding factors are pointing in different directions or too little information was available.
3. Unlikely: The number and/or specificity of the confounding factors suggest that a causal relationship is unlikely.

A total of three serious cases were received. In two of the serious cases it was evaluated that the experienced adverse events (AEs) were not related to the observed broken/cracked snaps on the cartridge holder while in the third limited and conflicting information was given. In the first case the reported AE did not match the medical consequence listed in the result code for the relevant fault and in the second case the patient consistently re-used needles and left them attached to the pen in between injections. Both needle re-use and leaving the needle attached to the pen are confounding factors which is more likely to cause the reported hyperglycaemic event. In the last case conflicting information related to the reported event and the patients' blood glucose measurements made the outcome of the assessment unknown.

In 18 non-serious cases the data available suggested that there was a causal relationship between the reported AEs and the observed fault(s) on the returned pen. In all of the 18 cases the patient experienced a hyperglycaemic event and in 5 of the 18 cases, the case was medical confirmed by a health care professional. Several factors which suggested a causal relationship or ruled out other potential explanations were present in the cases, please note that more than one factor can be present in one case.

<p>The HHE is a class II: A situation in which the use of or exposure to a violative product may cause temporary or medically reversible adverse health consequences, or where the probability of serious adverse health consequences is remote.</p>	
<p>Description and justification of the action (corrective / preventive)</p> <p>Corrective actions have been implemented in production 1st of Sep 2016 ensuring that all new NovoPen® Echo® and NovoPen® 5 devices have a cartridge holder without this potential issue.</p> <p>The consequence of a patient level recall will be a stock out situation for at least 4-6 weeks in certain markets which will put patients in a poor situation as it can have health consequences to switch to another injection device or another type of medication. To ensure patient safety it is decided to:</p> <p>1) provide affected patients with a cartridge holder of the original and more robust material, and which they can easily replace themselves. Novo Nordisk is establishing the logistics enabling courier service distribution of exchange cartridge holders to all affected patients globally. The cartridge holder is a part of the device and the patient is used to detach the cartridge holder, load a new insulin cartridge and attach the cartridge holder again.</p> <p>2) recall affected devices from wholesalers, pharmacies, hospitals and doctors.</p> <p>For clinical trial NN1218-4101, investigators and patients in the possession of affected NovoPen® Echo® devices will be contacted to have their devices replaced.</p>	
<p>Advice on actions to be taken by the distributor and the user</p> <p>Patients/user:</p> <ul style="list-style-type: none"> • Contact Novo Nordisk which will send a new, unaffected cartridge holder for the NovoPen® Echo® and/or NovoPen® 5 devices. The new, exchanged cartridge holder should be attached and used as stated in the Instruction For Use. <p>Doctors:</p> <p>Return affected devices to Novo Nordisk and advice patients to contact Novo Nordisk to have the cartridge holder replaced</p> <p>Wholesaler and pharmacies:</p> <ul style="list-style-type: none"> • Do not sell or hand out more of the affected devices. • Return affected devices to Novo Nordisk and order new devices from Novo Nordisk. 	
<p>Progress of FSCA , together with reconciliation data (Mandatory for a Final FSCA)</p>	
<p>Time schedule for the implementation of the different actions</p> <p>2017-07-05: Initiate Field Safety Corrective Actions outlined above in agreement with local health authorities</p>	
<p>Attached please find</p> <p><input checked="" type="checkbox"/> Field Safety Notice (FSN) in English</p> <p><input type="checkbox"/> FSN in national language</p> <p><input type="checkbox"/> Others (please specify)</p>	<p>FSN Status</p> <p><input type="radio"/> Draft FSN</p> <p><input checked="" type="radio"/> Final FSN</p>

The medical device has been distributed to the following countries:

within the EEA and Switzerland

- | | | | | | | | |
|--|--|--|--|--|--|--|--|
| <input checked="" type="checkbox"/> AT | <input checked="" type="checkbox"/> BE | <input checked="" type="checkbox"/> BG | <input checked="" type="checkbox"/> CH | <input type="checkbox"/> CY | <input checked="" type="checkbox"/> CZ | <input checked="" type="checkbox"/> DE | <input checked="" type="checkbox"/> DK |
| <input type="checkbox"/> EE | <input checked="" type="checkbox"/> ES | <input checked="" type="checkbox"/> FI | <input checked="" type="checkbox"/> FR | <input checked="" type="checkbox"/> GB | <input checked="" type="checkbox"/> GR | <input checked="" type="checkbox"/> HU | <input checked="" type="checkbox"/> IE |
| <input checked="" type="checkbox"/> IS | <input checked="" type="checkbox"/> IT | <input type="checkbox"/> LI | <input type="checkbox"/> LT | <input type="checkbox"/> LU | <input type="checkbox"/> LV | <input type="checkbox"/> MT | <input checked="" type="checkbox"/> NL |
| <input checked="" type="checkbox"/> NO | <input checked="" type="checkbox"/> PL | <input checked="" type="checkbox"/> PT | <input checked="" type="checkbox"/> RO | <input checked="" type="checkbox"/> SE | <input checked="" type="checkbox"/> SI | <input checked="" type="checkbox"/> SK | <input type="checkbox"/> TR |

Candidate Countries

- HR

- All EEA, candidate countries and Switzerland

Others:

Marketed products: Australia, Austria, Belgium, Bulgaria, Canada, Croatia, Czech Republic, Denmark, Finland, France, Georgia, G

8 Comments

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorised representative or the National Competent Authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.

Signature *2017-07-25 Eric Stephens*

I affirm that the information given above is correct
to the best of my knowledge

print

check

send XML-data by E-Mail

List of NovoPen® Echo® batches

Batch number as printed on the Pen

Country	Brand Name	Batch number
Australia	NovoPen® Echo®	EVG3310
Australia	NovoPen® Echo®	EVG4252
Australia	NovoPen® Echo®	EVG4253
Australia	NovoPen® Echo®	FVG8218
Australia	NovoPen® Echo®	FVG8412
Australia	NovoPen® Echo®	FVG8413
Australia	NovoPen® Echo®	FVG8415
Austria	NovoPen® Echo®	DUG2055
Austria	NovoPen® Echo®	DUG2058
Austria	NovoPen® Echo®	DVG1368
Austria	NovoPen® Echo®	DVG1564
Austria	NovoPen® Echo®	DVG1566
Austria	NovoPen® Echo®	DVG1715
Austria	NovoPen® Echo®	EVG5699
Austria	NovoPen® Echo®	EVG5946
Austria	NovoPen® Echo®	EVG6824
Austria	NovoPen® Echo®	FVG7337
Austria	NovoPen® Echo®	FVG8131
Austria	NovoPen® Echo®	FVG8132
Austria	NovoPen® Echo®	FVG8414
Austria	NovoPen® Echo®	FVG8995
Austria	NovoPen® Echo®	FVG8997
Belgium	NovoPen® Echo®	EVG3309
Belgium	NovoPen® Echo®	EVG3310
Belgium	NovoPen® Echo®	EVG4140
Belgium	NovoPen® Echo®	EVG5698
Belgium	NovoPen® Echo®	EVG5699
Belgium	NovoPen® Echo®	EVG6823
Belgium	NovoPen® Echo®	EVG6824
Belgium	NovoPen® Echo®	FVG7337
Bulgaria	NovoPen® Echo®	EVG2915
Bulgaria	NovoPen® Echo®	EVG5699
Bulgaria	NovoPen® Echo®	FVG8132
Bulgaria	NovoPen® Echo®	FVG8414
Canada	NovoPen® Echo®	DUG2373
Canada	NovoPen® Echo®	DVG1368
Canada	NovoPen® Echo®	EVG2915
Canada	NovoPen® Echo®	EVG3010
Canada	NovoPen® Echo®	EVG4252
Canada	NovoPen® Echo®	EVG5699
Canada	NovoPen® Echo®	EVG5946
Canada	NovoPen® Echo®	EVG5963
Canada	NovoPen® Echo®	FVG7337
Canada	NovoPen® Echo®	FVG7561
Canada	NovoPen® Echo®	FVG8131
Canada	NovoPen® Echo®	FVG8132

List of NovoPen® Echo® batches

Batch number as printed on the Pen

Canada	NovoPen® Echo®	FVG8212
Canada	NovoPen® Echo®	FVG8414
Canada	NovoPen® Echo®	FVG8994
Canada	NovoPen® Echo®	FVG8997
Croatia	NovoPen® Echo®	EVG3309
Croatia	NovoPen® Echo®	FVG7337
Croatia	NovoPen® Echo®	FVG7457
Croatia	NovoPen® Echo®	FVG8216
Czech Republic	NovoPen® Echo®	DVG1566
Czech Republic	NovoPen® Echo®	DVG1716
Czech Republic	NovoPen® Echo®	EVG2914
Czech Republic	NovoPen® Echo®	EVG2915
Czech Republic	NovoPen® Echo®	FVG8132
Czech Republic	NovoPen® Echo®	FVG8216
Czech Republic	NovoPen® Echo®	FVG8217
Czech Republic	NovoPen® Echo®	FVG8218
Czech Republic	NovoPen® Echo®	FVG8412
Denmark/Iceland	NovoPen® Echo®	DUG1709
Denmark/Iceland	NovoPen® Echo®	DUG1778
Denmark/Iceland	NovoPen® Echo®	DUG2116
Denmark/Iceland	NovoPen® Echo®	DUG2130
Denmark/Iceland	NovoPen® Echo®	EVG2298
Denmark/Iceland	NovoPen® Echo®	EVG2300
Denmark/Iceland	NovoPen® Echo®	EVG4252
Denmark/Iceland	NovoPen® Echo®	EVG4253
Denmark/Iceland	NovoPen® Echo®	EVG5699
Denmark/Iceland	NovoPen® Echo®	EVG5946
Denmark/Iceland	NovoPen® Echo®	EVG6823
Denmark/Iceland	NovoPen® Echo®	EVG6824
Denmark/Iceland	NovoPen® Echo®	FVG8131
Denmark/Iceland	NovoPen® Echo®	FVG8132
Denmark/Iceland	NovoPen® Echo®	FVG8413
Denmark/Iceland	NovoPen® Echo®	FVG8414
Finland	NovoPen® Echo®	DUG2116
Finland	NovoPen® Echo®	DUG2130
Finland	NovoPen® Echo®	DUG2219
Finland	NovoPen® Echo®	DUG2343
Finland	NovoPen® Echo®	DUG2431
Finland	NovoPen® Echo®	DUG2436
Finland	NovoPen® Echo®	DVG1564
Finland	NovoPen® Echo®	DVG1565
Finland	NovoPen® Echo®	DVG2297
Finland	NovoPen® Echo®	EVG2299
Finland	NovoPen® Echo®	EVG3309
Finland	NovoPen® Echo®	EVG3310
Finland	NovoPen® Echo®	EVG4252
Finland	NovoPen® Echo®	EVG5699

List of NovoPen® Echo® batches

Batch number as printed on the Pen

Finland	NovoPen® Echo®	EVG6823
Finland	NovoPen® Echo®	EVG6824
Finland	NovoPen® Echo®	FVG7561
Finland	NovoPen® Echo®	FVG7571
Finland	NovoPen® Echo®	FVG8414
Finland	NovoPen® Echo®	FVG8997
France	NovoPen® Echo®	DUG2431
France	NovoPen® Echo®	DUG2436
France	NovoPen® Echo®	DVG2297
France	NovoPen® Echo®	EVG2299
France	NovoPen® Echo®	EVG2915
France	NovoPen® Echo®	EVG3010
France	NovoPen® Echo®	EVG3309
France	NovoPen® Echo®	EVG5946
France	NovoPen® Echo®	EVG6824
France	NovoPen® Echo®	FVG7457
France	NovoPen® Echo®	FVG7561
France	NovoPen® Echo®	FVG7571
France	NovoPen® Echo®	FVG8132
France	NovoPen® Echo®	FVG8994
France	NovoPen® Echo®	FVG8995
France	NovoPen® Echo®	FVG8997
Georgia	NovoPen® Echo®	EVG2915
Georgia	NovoPen® Echo®	EVG3309
Georgia	NovoPen® Echo®	FVG8132
Georgia	NovoPen® Echo®	FVG8217
Germany	NovoPen® Echo®	DUG1776
Germany	NovoPen® Echo®	DUG2052
Germany	NovoPen® Echo®	DUG2053
Germany	NovoPen® Echo®	DUG2056
Germany	NovoPen® Echo®	DUG2057
Germany	NovoPen® Echo®	DUG2058
Germany	NovoPen® Echo®	DUG2116
Germany	NovoPen® Echo®	DUG2130
Germany	NovoPen® Echo®	DUG2344
Germany	NovoPen® Echo®	DUG2431
Germany	NovoPen® Echo®	DUG2432
Germany	NovoPen® Echo®	DUG2436
Germany	NovoPen® Echo®	DVG1368
Germany	NovoPen® Echo®	DVG1564
Germany	NovoPen® Echo®	DVG1564
Germany	NovoPen® Echo®	DVG1716
Germany	NovoPen® Echo®	EVG2300
Germany	NovoPen® Echo®	EVG2908
Germany	NovoPen® Echo®	EVG2915
Germany	NovoPen® Echo®	EVG3309
Germany	NovoPen® Echo®	EVG3640

List of NovoPen® Echo® batches

Batch number as printed on the Pen

Germany	NovoPen® Echo®	EVG3641
Germany	NovoPen® Echo®	EVG4140
Germany	NovoPen® Echo®	EVG4252
Germany	NovoPen® Echo®	EVG5698
Germany	NovoPen® Echo®	EVG5716
Germany	NovoPen® Echo®	EVG5946
Germany	NovoPen® Echo®	EVG5964
Germany	NovoPen® Echo®	FVG7149
Germany	NovoPen® Echo®	FVG7457
Germany	NovoPen® Echo®	FVG7458
Germany	NovoPen® Echo®	FVG8134
Germany	NovoPen® Echo®	FVG8135
Germany	NovoPen® Echo®	FVG8412
Germany	NovoPen® Echo®	FVG8413
Greece	NovoPen® Echo®	EVG5946
Greece	NovoPen® Echo®	FVG8413
Hong Kong	NovoPen® Echo®	EVG6823
Hong Kong	NovoPen® Echo®	FVG8413
Hungary	NovoPen® Echo®	EVG6823
Hungary	NovoPen® Echo®	EVG6824
Italy	NovoPen® Echo®	DUG2374
Italy	NovoPen® Echo®	DVG2297
Italy	NovoPen® Echo®	EVG2914
Italy	NovoPen® Echo®	FVG8131
Italy	NovoPen® Echo®	FVG8216
Italy	NovoPen® Echo®	FVG8997
Kyrgyzstan	NovoPen® Echo®	DVG1368
Kyrgyzstan	NovoPen® Echo®	DVG1564
Kyrgyzstan	NovoPen® Echo®	EVG5946
Netherlands	NovoPen® Echo®	DUG2343
Netherlands	NovoPen® Echo®	DUG2344
Netherlands	NovoPen® Echo®	DVG1566
Netherlands	NovoPen® Echo®	DVG1715
Netherlands	NovoPen® Echo®	EVG5963
Netherlands	NovoPen® Echo®	EVG6823
Netherlands	NovoPen® Echo®	EVG6824
Netherlands	NovoPen® Echo®	FVG8413
Netherlands	NovoPen® Echo®	FVG8414
New Zealand	NovoPen® Echo®	FVG8412
New Zealand	NovoPen® Echo®	FVG8415
Norway	NovoPen® Echo®	DUG1709
Norway	NovoPen® Echo®	DUG1778
Norway	NovoPen® Echo®	DVG1368
Norway	NovoPen® Echo®	DVG1564
Norway	NovoPen® Echo®	EVG2915
Norway	NovoPen® Echo®	EVG3011
Norway	NovoPen® Echo®	EVG3310

List of NovoPen® Echo® batches

Batch number as printed on the Pen

Norway	NovoPen® Echo®	EVG6823
Norway	NovoPen® Echo®	EVG6824
Norway	NovoPen® Echo®	FVG7561
Norway	NovoPen® Echo®	FVG8413
Norway	NovoPen® Echo®	FVG8414
Norway	NovoPen® Echo®	FVG8415
Norway	NovoPen® Echo®	FVG8995
Poland	NovoPen® Echo®	EVG2299
Poland	NovoPen® Echo®	EVG2914
Poland	NovoPen® Echo®	EVG3011
Poland	NovoPen® Echo®	EVG3310
Poland	NovoPen® Echo®	EVG3545
Poland	NovoPen® Echo®	EVG4140
Poland	NovoPen® Echo®	FVG7364
Poland	NovoPen® Echo®	FVG7457
Poland	NovoPen® Echo®	FVG7561
Poland	NovoPen® Echo®	FVG7571
Poland	NovoPen® Echo®	FVG8131
Poland	NovoPen® Echo®	FVG8132
Poland	NovoPen® Echo®	FVG8413
Poland	NovoPen® Echo®	FVG8414
Portugal	NovoPen® Echo®	EVG3010
Portugal	NovoPen® Echo®	FVG7561
Portugal	NovoPen® Echo®	FVG8132
Portugal	NovoPen® Echo®	FVG8414
Romania	NovoPen® Echo®	EVG5699
Romania	NovoPen® Echo®	EVG5963
Romania	NovoPen® Echo®	FVG8218
Russia	NovoPen® Echo®	FVG8131
Russia	NovoPen® Echo®	FVG8216
Russia	NovoPen® Echo®	FVG8413
Russia	NovoPen® Echo®	FVG8996
Serbia	NovoPen® Echo®	EVG2299
Serbia	NovoPen® Echo®	EVG5699
Singapore	NovoPen® Echo®	FVG8413
Slovakia	NovoPen® Echo®	EVG4252
Slovakia	NovoPen® Echo®	EVG4253
Slovenia	NovoPen® Echo®	EVG2298
Slovenia	NovoPen® Echo®	FVG8997
South Africa	NovoPen® Echo®	DVG1368
South Africa	NovoPen® Echo®	DVG1565
Spain	NovoPen® Echo®	DVG1564
Spain	NovoPen® Echo®	EVG2298
Spain	NovoPen® Echo®	EVG2914
Spain	NovoPen® Echo®	EVG3011
Spain	NovoPen® Echo®	EVG4140
Spain	NovoPen® Echo®	FVG7571

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