Johnson & Johnson LLC

1-V, Pavla Tychyny avenue, Kiev, 02152
Tel.: +38 (044) 498 08 88 Fax: +38 (044) 498 73 92

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(MEDDEV 2.12/1 rev 7)

Version 2.7en
2012-12-03

1 Administrative information
To which NCA(s) is this report being sent? Bundesinstitut für Arzneimittel und Medizinprodukte Dr. Ekkehard Stöblein Kurt Georg Kiesinger Allee 3 D - 53175 Bonn
Type of report <input type="radio"/> Initial report <input type="radio"/> Follow-up report <input checked="" type="radio"/> Final report
Date of this report 2019-06-27
Reference number assigned by the manufacturer Harmonic ACE+ (HAR23 and HAR36) - VOLUNTARY RECALL
FSCA reference number assigned by NCA DE-BfArM-2018-05-17-2476
Incidence reference number assigned by NCA 05454/18
Name of the co-ordinating NCA Competent Authority (if applicable) BfArM

2 Information on submitter of the report
Status of submitter <input type="radio"/> Manufacturer <input checked="" type="radio"/> Authorised Representative within EEA and Switzerland <input type="radio"/> Others: (identify the role)

3 Manufacturer information	<input type="button" value="new"/>
Name Ethicon Endo Surgery , LLC	
Contact Name Nathan Anderson	
Address 475 C Street, Los Frailes Industrial Park, Suite 401	
Postcode 00969	City Guaynabo, Puerto Rico
Phone 0017872726144	Fax 0017872723003
E-mail NAnders1@its.jnj.com	Country US - USA

4 Authorised Representative Information

new

Name Ethicon Endo Surgery (Europe) GmbH	
Contact Name Jörg Hansel	
Address Hummelsbüttler Steindamm 71	
Postcode 22851	City Norderstedt
Phone 00494052974232	Fax
E-mail jhansel@its.jnj.com	Country DE - Germany

5 National contact point information

new

National contact point name see attached contact list	
Name of the contact person see attached contact list	
Address see attached contact list	
Postcode see attached contact list	City see attached contact list
Phone see attached contact list	Fax
E-mail see attached contact list	Country DE - Germany

6 Medical device information

new

Class

- AIMD Active implants
 MDD Class III
 MDD Class IIb
 MDD Class IIa
 MDD Class I
- IVD Annex II List A
 IVD Annex II List B
 IVD Devices for self-testing
 IVD General

Nomenclature system (preferable GMDN)

GMDN

Nomenclature code

44755

Nomenclature text

Ultrasonic Surgical System Handpiece Tip

Commercial name/ brand name / make

See product details in attached list

Model number

N/A

Catalogue number

See product details in attached list

Serial number(s)

N/A

Lot/batch number(s)

See product details in attached list

Device Mfr Date

2016-10-31

Expiry date

2022-11-23

Notified Body (NB) ID-number

0123

Accessories / associated devices (if applicable)

N/A

Software version number (if applicable)

N/A

7 Description of the FSCA

Background information and reason for the FSCA

Recently, Ethicon has become aware of reports from customers related to continuouse/inadvertent activation during analysis with HARMONIC ACE® LAPAROSCOPIC 5MM DIAMETER SHEARS + ADAPTIVE TISSUE TECHNOLOGY (HAR23 and HAR36).

Based on preliminary investigation, the potential root cause for the inverted dome is believed to be damage to the flex circuit assembly caused at the Plexus (supplier) manufacturing/inspection process or due to dimensional design differences, causing assembly complications at the Juarez manufacturing plant. Ethicon had placed stop shipments on the subject lots of HARMONIC ACE® LAPAROSCOPIC 5MM DIAMETER SHEARS + ADAPTIVE TISSUE TECHNOLOGY and instituted this medical device recall.

Description and justification of the action (corrective / preventive)

Ethicon decided to initiate a voluntary medical device recall of certain lots of HARMONIC ACE® LAPAROSCOPIC 5MM DIAMETER SHEARS + ADAPTIVE TISSUE TECHNOLOGY (HAR23 and HAR36) distributed to the field.

Based on preliminary investigation, the potential root cause for the inverted dome is believed to be damage to the flex circuit assembly caused at the Plexus (supplier) manufacturing/inspection process or due to dimensional design differences, causing assembly complications at the Juarez manufacturing plant. Ethicon had placed stop shipments on the subject lots of HARMONIC ACE® LAPAROSCOPIC 5MM DIAMETER SHEARS + ADAPTIVE TISSUE TECHNOLOGY and instituted this medical device recall.

Ethicon had stopped manufacturing HARMONIC ACE® devices with Plexus components as of December 12, 2017.

A CAPA will be originated to further investigate the issue, root cause, and implement corrective and/or preventive actions.

Advice on actions to be taken by the distributor and the user

See attached FSN

Progress of FSCA , together with reconciliation data (Mandatory for a Final FSCA)

- Initial FSCA sent on April 26th, 2018.
- First Follow Up FSCA sent on October 5th, 2018.
- Second Follow Up FSCA sent on April 01th, 2019
- Final FSCA sent on June, 27th, 2019
- CAPA-008265: All CAPA Action completed
- All customers within the scope of the voluntary recall were identified and notified.
- All available products were returned to Ethicon

Time schedule for the implementation of the different actions

- Initial FSCA sent on April 26th, 2018.
- First Follow Up FSCA sent on October 5th, 2018.
- Second Follow Up FSCA sent on April 01th, 2019
- Final FSCA sent on June, 27th, 2019

Attached please find

- Field Safety Notice (FSN) in English
- FSN in national language
- Others (please specify)

see Email

FSN Status

- Draft FSN
- Final FSN

The medical device has been distributed to the following countries:

within the EEA and Switzerland

- | | | | | | | | |
|--|--|--|--|--|--|--|--|
| <input type="checkbox"/> AT | <input checked="" type="checkbox"/> BE | <input checked="" type="checkbox"/> BG | <input checked="" type="checkbox"/> CH | <input type="checkbox"/> CY | <input checked="" type="checkbox"/> CZ | <input checked="" type="checkbox"/> DE | <input checked="" type="checkbox"/> DK |
| <input checked="" type="checkbox"/> EE | <input checked="" type="checkbox"/> ES | <input checked="" type="checkbox"/> FI | <input checked="" type="checkbox"/> FR | <input checked="" type="checkbox"/> GB | <input checked="" type="checkbox"/> GR | <input checked="" type="checkbox"/> HU | <input checked="" type="checkbox"/> IE |
| <input checked="" type="checkbox"/> IS | <input checked="" type="checkbox"/> IT | <input type="checkbox"/> LI | <input checked="" type="checkbox"/> LT | <input checked="" type="checkbox"/> LU | <input checked="" type="checkbox"/> LV | <input checked="" type="checkbox"/> MT | <input checked="" type="checkbox"/> NL |
| <input checked="" type="checkbox"/> NO | <input checked="" type="checkbox"/> PL | <input checked="" type="checkbox"/> PT | <input checked="" type="checkbox"/> RO | <input checked="" type="checkbox"/> SE | <input checked="" type="checkbox"/> SI | <input checked="" type="checkbox"/> SK | <input checked="" type="checkbox"/> TR |

Candidate Countries

- HR

- All EEA, candidate countries and Switzerland

Others:

Algeria, Angola, Argentina, Armenia, Australia, Bolivia, Bosnia and Herzegovina, Brazil, Canada, Chile, China, Colombia, Costa Rica+

8 Comments

N/A

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorised representative or the National Competent Authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.

Signature

Milan
Miladinovic

I affirm that the information given above is correct to the best of my knowledge

print

check

send XML-data by E-Mail