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URGENT – Field Safety Notice

Commercial name of the affected product	GMK Sphere Flex Inserts
FSCA-identifier	FSCA 2024-01 – 25/09/2024
Type of action	Recall of two specific lots
SRN – single registration number	CH-MF-000033473

Affected Products

Description	UDI	Reference	Lot
GMK-SPHERE Tibial Insert - Flex S6L - 12 mm	07630030826979	02.12.0612FL	2337090
GMK-SPHERE Tibial Insert - Flex S6L - 11 mm	07630030862687	02.12.0611FL	2337089

1. Description of advisory:

Medacta International is conducting a removal of two lots from the market for the following reason:

During a knee surgery, a tibial insert size 6, 11mm, reference 02.12.0611FL, lot 2337089 was found inside the box of a tibial insert size 6, 12mm, labelled as reference 02.12.0612FL, lot 2337090.

Upon internal investigation, a mix-up of the two lots has been confirmed.

The use of an insert that is 1mm thicker or thinner than the trial (11 mm instead of 12 mm and vice versa) could marginally affect joint stability, soft tissue tension, and potentially range of motion, potentially resulting in some discomfort for the patient. Insert thickness can be increased by an increment of 1 or 2mm based on the state-of-the-art: therefore, although an unintended offset of 1mm may be perceivable, it is not likely to impair the outcome. Checking the stability and soft tissue tension are part of the standard procedure prior to completing a TKA surgery, therefore a situation that could lead to an unsatisfactory condition is likely to be detected by the operating surgeon. Based on these factors, no pre-emptive interventions are recommended, however monitoring of the joint stability and subjective feeling should be evaluated during the post-operative patient follow-up visits, as per the standard of care.

2. Action to take:

- Return all units of the lot numbers listed to Medacta International SA
- Forward this notice to all people and organizations potentially affected.

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Strada Regina
CH-6874 Castel San Pietro-Switzerland

Phone +41 91 696 60 60
Fax +41 91 696 60 66

Medacta confirms that all relevant Competent Authorities and Regulatory Agencies have been advised of this Field Safety Notice
This may be translated, however the legally binding version is the original English issued by Medacta International SA

3. Contact person and return address:

Stefano Baj, Regulatory Affairs Director
Medacta International SA, Strada Regina, CH-6847 Castel San Pietro, Switzerland
Tel: +41 (0)91 696 60 60, Fax: +41 (0)91 696 60 66

We apologize for this inconvenience and confirm our commitment to patient safety and surgeon support.

Sincerely,

Stefano Baj
Regulatory Affairs Director

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